



Suma Med Rinse

Rinse aid for (invasive and non-invasive) medical devices

Description

Suma Med Rinse is a concentrated low foam rinse additive for (invasive and non-invasive) medical devices cleaning in machine warewashing.

Key properties

Suma Med Rinse is a concentrated neutral rinse additive. It has been specially formulated for use in a wide range of warewashing machines. The special blend of non-ionic surfactants is designed to lower the surface tension of the final rinsing water. This prevents droplet formation, and because the water sheets off the surface in a thin layer, drying times are also reduced. This special blend also prevents heavy foaming in the washer disinfectant machine.

Benefits

- Ensures quick drying
- Gives streak and spot-free results
- Its concentrated formulation gives excellent cost in use
- Low foaming action for a more efficient operation

Use instructions

Suma Med Rinse is recommended for use in all types of automatically dosed commercial ware washing machines. The product is directly injected into the final rinse water, using automatic dosing equipment. Dry medical device before use. Additionally, medical devices would be sterilized further if they are invasive. Safety glasses are not normally required. However, their use is recommended in those cases where splashes may occur when handling the product

Dosage*

Use at a concentration of 0.2-0.5ml/L at temperatures of 70-80°C in soft/medium water conditions. Actual dosage levels are determined according to local site conditions.

**This dosage is according to optimal conditions, recommendations may vary, please consult with your Diversey representative for directions.*

Application

Suma Med Rinse is intended to rinse non-invasive and invasive medical devices after being cleaned before medical treatment. It also makes the drying process more efficient. The product is intended for single-use. The product is used in the final rinse water of a washing machine of reusable medical devices. After drying, there are no expected residues that come into contact with the skin of the patient.

This product is an accessory to the following medical devices:


- Suma Med Super LpH
- Suma Med Super Plus
- Suma Med Alu L10
- Suma Med Enzyme

Intended Users

Authorized trained personnel at a healthcare facility. For professional use only.



 Diversey Europe
Operations B.V.
Maarssebroeksedijk 2,
3542 DN Utrecht,
The Netherlands





Suma Med Rinse

Technical data

Appearance: Clear, colorless liquid
pH value (neat): ≈ 7
pH value (in use): ≈ 7 (0.05 %)
Relative density (g/cm³; 20°C): ≈ 1.02 (20 °C)

The above data is typical of normal production and should not be taken as a specification.

Safe handling and storage information

Full guidance on the handling and disposal of this product is provided in a separate Safety Data Sheet; sds.diversey.com. Handle following good industrial hygiene and safety practices. Keep away from food, drink, and animal feedstuff. Do not mix with other products unless advised by Diversey. Store in original closed containers away from extremes of temperature. The product is not intended to be re-used. Only for professional users/specialists.

Warning

Storage temperatures 0-40 °C.

In case of any serious incident occurs in relation to the use of Suma Med Rinse, the user should report it to local Diversey office and to the competent authority of the Member State in which the user and/or patient is established. Do not re-use. The product is not intended for use on the human skin. Not intended for soaking. Suma Med Rinse is not intended to be used on porous. Do not use on Medical Devices made of Polyethylene Terephthalate (PET). Without proper sterilization, the product is not intended for use on invasive medical devices that penetrate tissue or invade the bloodstream nor come into direct contact with intact mucous membranes. Uses other than those identified are not recommended. Wash hands before breaks and at the end of the workday. Wash face, hands, and any exposed skin thoroughly after handling. Take off immediately all contaminated clothing. Use personal protective equipment as required. Use only with adequate ventilation. May produce an allergic reaction.

Safety data sheet available on request. Do not use Suma Med Rinse on plastic medical devices.

Personal Protective Equipment

Recommended safety measures for handling the undiluted product: Covering activities such as filling and transfer of product to application equipment, flasks, or buckets.

| Personal Protective Equipment | |
|-------------------------------|---|
| Eye/face protection | Safety glasses are not normally required. However, their use is recommended in those cases where splashes may occur when handling the product (EN 166). |

Recommended safety measures for handling the diluted product: No special requirements under normal use conditions.

Disposal

The concentrated contents or contaminated packaging should be disposed of by a certified handler or according to the site permit. The release of waste to sewers is discouraged. The cleaned packaging material is suitable for energy recovery or recycling in line with local legislation.

European Waste Catalogue: 20 01 30 - detergents other than those mentioned in 20 01 29.

Empty packaging Recommendation: Dispose of observing national or local regulations.

Suitable cleaning agents: Water, if necessary with a cleaning agent.

www.diversey.com

© 2023 Diversey, Inc. All Rights Reserved. 07/04/2023 en (P01157)